#### (19) World Intellectual Property Organization International Bureau



### 

(43) International Publication Date 11 January 2001 (11.01.2001)

PCT

# (10) International Publication Number

WO 01/01925 A1

- (51) International Patent Classification7: C03C 10/00
- A61K 6/02,
- (74) Agent: REIMER, Leah, M.; Cantor Colburn LLP, 55 Griffin Road South, Bloomfield, CT 06002 (US).

- (21) International Application Number:
- PCT/US00/40296
- (22) International Filing Date:
- 30 June 2000 (30.06.2000)
- (25) Filing Language:

English

(26) Publication Language:

English

(30) Priority Data: 60/142,204

2 July 1999 (02.07.1999)

- (71) Applicant: JENERIC/PENTRON INCORPORATED [US/US]; 53 North Plains Industrial Road, Wallingford, CT 06492 (US).
- (72) Inventors: BRODKIN, Dmitri; 75 Herbert Terrace, West Orange, NJ 07052 (US). PANZERA, Carlino; 7 Huntsman Lane, Bellemead, NJ 07052 (US). PANZERA, Paul; 1 Middletown Drive, Mt. Holly, NJ 08060 (US).

(81) Designated State (national): JP.

(84) Designated States (regional): European patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE).

#### Published:

- With international search report.
- Before the expiration of the time limit for amending the claims and to be republished in the event of receipt of amendments

For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.



(54) Title: DENTAL PORCELAINS

(57) Abstract: Opaque porcelains for use with metal cores in the manufacture of PFM restorations. The porcelains exhibit a coefficient of thermal expansion (CTE) substantially equal to or slightly above the CTE of the metal to which it is applied. The porcelains exhibit a CTE equal to or up to about 1.5 x 10-6/°C higher than the dental alloys to which they are applied as the opaque. The porcelains are fabricated from a mixture of two frit compositions. A high expansion, leucite containing frit is combined with a low melting glass frit to provide a porcelain having an expansion in the range of 16.9 to about 18 x 10% C in the temperature range of 25° - 500°C.

WO 01/01925 PCT/US00/40296

#### **DENTAL PORCELAINS**

#### Field of the Invention

This invention relates to low fusing high-expansion dental porcelain especially useful for the fabrication of both all-ceramic and porcelain-fused-to-metal (PFM) restorations.

### Background of the Invention

5

10

15

Porcelains are typically designed to be used in the manufacture of either all-ceramic dental restorations or in PFM restorations, but are not normally functional with both types of restorations due to the differences in properties of ceramics and metals. One such porcelain, OPC® Low Wear™ porcelain, available from Jeneric/Pentron Inc., Wallingford, CT and covered in copending, commonly assigned patent application No. 09/133,582 filed August 13, 1998 which is hereby incorporated by reference, was initially designed to be used as overlay for pressed all-ceramic restorations as well as for the fabrication of porcelain jacket crowns and veneers. However, OPC®Low Wear™ porcelain is not currently used for PFM restorations despite its wear resistance, forgiveness to natural dentition and strength being superior to those of conventional PFM porcelains as shown in the Table 1 below:

15

Table 1

Property	OPC <sup>®</sup> Low Wear <sup>™</sup> Porcelain	Conventional Porcelain for PFM
Leucite average grain Size, µm	About 2-3	About 5-8
Leucite volume fraction, %	35-40	20-25
Enamel wear*, x10 <sup>-2</sup> mm <sup>2</sup>	7.69 <u>+</u> 3.20	18.23 <u>+</u> 5.20
Wear of ceramics*, x10 <sup>-2</sup> mm <sup>3</sup>	0.16 <u>+</u> 0.04	0.49 <u>+</u> 0.11
CTE, 10 <sup>-6</sup> /°C, 25°C-500°C	About 17	About 13

The major obstacle preventing use of the OPC® Low Wear<sup>TM</sup> porcelain in

5 PFM restorations is the absence of an opaque/alloy combination compatible with this porcelain having relatively high expansion of about 17x10<sup>-6</sup>/°C (25°C-500°C).

There exists a Golden Gate System<sup>TM</sup> for PFM restorations available through Degussa<sup>TM</sup> (Dental Division, South Plainfield, NJ) which combines Duceragold<sup>TM</sup> porcelain and Degunorn<sup>TM</sup>, type IV crown and bridge alloy (CTE=16.4x10<sup>-6</sup>/°C, 25°C-500°C). This system requires rather tedious multistep alloy preparation procedures including a necessary wash bake step prior to application of the opaque; and an excessively long (16-20 min) and complex first dentine bake to assure proper bonding and compatibility of the Duceragold<sup>TM</sup> porcelain to the Degunorm<sup>TM</sup> alloy. In particular, the cooling segment (3-4 min between 720°C and 680°C) in the first dentine bake is required by the manufacturer to grow additional leucite and may be an indication of instability of leucite in this porcelain. The following Table 2 below sets forth the various properties of the Duceragold<sup>TM</sup> porcelain.

Table 2

	Duceragold™
Firing Temperature, °C	770-790
Glass Transition	490
Temperature, °C	
Softening Temperature, °C	595
CTE 25° -500°C, 10°6/°C	15.8
Recommended alloy	Degunorm
Alloy CTE 25° -500°C, 10 <sup>-6</sup> /°C	16.4

10

15

20

30

U.S. Patent No. 5,453,290 to Van der Zel is directed to a dental porcelain for use with a dental alloy. The porcelain described therein must be fabricated from three frits, making it more difficult and costly to control the expansion and the glass transition temperature of the final product. Moreover, the CTE of the porcelain must be below the CTE of the alloy by 0.5-1.5 limiting the components to be used together. There is a need to provide a porcelain-fused-to-metal system for dental restorations having simple manufacturing procedures. It is desirable to provide a porcelain that is compatible with alloys of relatively high expansion. It is desirable to provide a two-frit porcelain for use in PFM restorations.

### Summary of the Invention

These and other objects and advantages are accomplished by opaque porcelains for use with metal cores in the manufacture of PFM restorations. The porcelains exhibit a coefficient of thermal expansion (CTE) substantially equal to or slightly above the CTE of the metal to which it is applied. In a preferred embodiment, the porcelains exhibit a CTE equal to or up to about  $1.0 \times 10^{-6}$ /°C higher than the dental alloys to which they are applied as the opaque. The porcelains are fabricated from a mixture of two frit compositions. A high expansion, leucite containing frit is combined with a low fusing glass frit to provide a porcelain having an expansion in the range of 16.9 to about  $18.5 \times 10^{-6}$ /°C in the temperature range of  $25^{\circ}$  -  $500^{\circ}$ C. By combining two frits, the expansion and fusing temperature can be controlled to the values stated above.

### 25 Description of the Invention

The invention relates to a porcelain material for use in all-ceramic restorations and PFM restorations. The porcelains exhibit a coefficient of thermal expansion (CTE) substantially equal to or slightly above the CTE of the metal to which it is applied. Preferably, the porcelains exhibit a CTE equal to or up to about 1.5 x 10<sup>-6</sup>/°C higher than the dental alloys to which they are applied as the opaque, and more

15

preferably equal to or up to about  $1.0 \times 10^{-6}$  C. The porcelains are compatible with metals having coefficients of thermal expansion (CTE) in the range of from about 15.5 to about 17 x 10-6/°C in the temperature range from 20° to 500 °C. The porcelains exhibit CTEs in the range of about 16.9 to about 18.5 in the temperature range from 20° to 500 °C, and preferably in the range of about 17 to about 17.5 x 10° <sup>6</sup>/°C in the temperature range from 20° to 500 °C. The porcelains are fabricated from a mixture of two frit compositions. A high expansion, leucite containing frit is combined with a low fusing glass frit to provide a porcelain having an expansion in the range of 16.9 to about 18.5  $\times$  10<sup>-6</sup>/°C in the temperature range of 25° - 500°C. It is essential to this invention that both the high expansion and the low fusing components of the two-frit mixture exhibit a low glass transition temperature (GTT). It is extremely important that the low fusing component used in the opaque formulation exhibits a GTT lower than about 415°C. By combining two frits, the expansion and firing temperature can be controlled to the values stated above. Opaque porcelains herein having pigments exhibit a coefficient of thermal expansion of average value of about 17 x 10<sup>-6</sup>/°C in the temperature range of 25° - 500°C. Opaque porcelains without pigments, i.e., white opaques, exhibit coefficients of thermal expansion in the higher end of the range, such as about 17.5 to about 18.5 x  $10^{-6}$ /°C in the temperature range of 25° - 500°C.

Table 3 below shows the compositional ranges of the porcelains for use in the invention.

1	`a	b	l	е	3
-	•	·	٠,	•	~

Table 5					
Body and	Opaque				
Incisal	Porcelain				
Porcelain					
59-65	59-65				
X	X				
10-15	10-15				
X	X				
0.5-2	0.5-2				
X	X				
X	X				
1.5-3	1.5-3				
15-17	15-17				
4-6	4-6				
X	X				
X	X				
X	х				
0.4-1	0.4-1				
	X				
0-1	13-20				
0-5	2-13				
	Body and Incisal Porcelain 59-65 X 10-15 X 0.5-2 X X 1.5-3 15-17 4-6 X X X 0.4-1 0-1				

<sup>\*</sup>Opacifiers, TiO<sub>2</sub>, ZrO<sub>2</sub>, ZrSiO<sub>4</sub>, ZnO, CeO<sub>2</sub>, or Ta<sub>2</sub>O<sub>5</sub>, are admixed as fine powder to a mixture of two frits. The resulting composition is referred to below as White Porcelain (opaque, body or incisal)

As set forth in Table 3 above, Li<sub>2</sub>O, present in an amount of from about 1.5% to about 3%, and F, present in an amount of 0.4%-1%, are instrumental in providing a low glass transition temperature. The presence of Li<sub>2</sub>O and F also assist as well in increasing the coefficient of thermal expansion and decreasing the maturing (firing) temperature. The high expansion, leucite containing component of the two-frit mixture has a reasonably low glass transition temperature as well. This is achieved by maintaining a reasonably low molar ratio of Al<sub>2</sub>O<sub>3</sub> to the sum of alkali and alkaline

<sup>\*\*</sup>Pigments are admixed as fine powder to a White Porcelain powder. The resulting powder composition is referred to below as Shaded Porcelain (opaque, body or incisal) X signifies non-essential components

15

earth oxides (R2O+RO). Normally, these compositions are extremely unstable and reactive as well as prone to sanidine precipitation in the temperature range of 650°C-950°C. It was surprisingly found that certain compositions with specific combinations of K<sub>2</sub>O, Na<sub>2</sub>O and Li<sub>2</sub>O are very stable. In addition, it was found that the molar ratio of 5 Al<sub>2</sub>O<sub>3</sub>/K<sub>2</sub>O should be within the range of 0.73-0.95 to assure both the required thermal stability and low glass transition temperature of the high expansion component of the porcelain. The low glass transition temperature provides a porcelain having good compatability with alloys having CTEs in the range of about 15.5 to about 17 x 10-6/°C in the temperature range of 25  $^{\circ}$  – 500  $^{\circ}$  C. Table 4 below sets forth the properties of the porcelain compositions.

Table 4

PORCELAIN	Body & Incisal	Opaque
Firing Temperature, °C	855-870	855-890
Glass Transition	420-430	430-440
Temperature, °C		
Softening Temperature, °C	520	
CTE 25° - 400°C, 10 <sup>-6</sup> /°C	15.2	
CTE 25° - GTT, 10°/°C	15.8 (GTT=430 °C)	<del></del>
CTE 25° -500°C, 10°/°C	17.2	17.0

Table 5 sets forth compatible alloys for use with the porcelains.

Table 5

1 4010 5				
Alloy	CTE (25°C-500°C)	Application		
Bio-75G	15.5±0.2	Fast Cool for single units		
		only		
GoldCore 75	16.4±0.2	For single units and bridges		
GoldCore 55	17.0±0.2	For single units and bridges		
JewelCast	17.0±0.2	For single units and bridges		

Table 6 below shows compositional examples of body (incisal) and opaque porcelains.

Table 6

			l able 6			
	Example 1	Comparative	Example 2	Example 3	Example 4	Example 5
	Body/Incisal	Example 1	White	White	White	White
	Porcelain	•	Opaque	Opaque	Opaque	Opaque
			Porcelain	Porcelain	Porcelain	Porcelain
			for Light*	for Light*	for Dark**	for Dark**
	!		Shades	Shades	Shades	Shades
Two-frit Mixt	ure Compos	tion:				
SiO2	61.9	61.5	60.0	58.5	60.1	58.6
B2O3	0.5	0.5	0.6	0.6	0.6	0.6
Al2O3	11.6	15.4	13.6	13.6	13.7	13.6
ZnO	0.0	0.0	0.0	2.0	0.0	2.0
CaO	1.7	0.6	1.2	1.2	1.2	1.2
MgO	0.8	0.2	0.2	0.1	0.2	0.1
BaO	0.0	0.0	0.4	0.4	0.4	0.4
Li2O	2.5	2.3	2.3	2.3	2.3	2.3
K20	15.7	12.9	15.1	15.0	15.1	15.0
Ta2O5	0.0		0.0	0.0	0.0	0.0
Na2O	4.8	6.0	5.5	5.4	5.4	5.3
TiO2	0.0	0.0	0.0	0.0	0.0	0.0
ZrO2	0.0	0.0	0.0	0.0	0.0	0.0
CeO2	0.0	0.6	0.3	0.3	0.3	0.3
F	0.6		0.7	0.6	0.7	0.6
Mixed-in Op		I				
ZrO2	T		17	17	15	15
ZrSiO4						
TiO2	† · · · · · · · · · · · · · · · · · · ·					
SnO2						
ZnO		1				
CeO2						
Ta2O5	<del> </del>					
Firing	857		871	871	871	871
temperature						
CTE (25°C-	17.2±0.3		17.8	17.8	17.8	17.8
500°C)						

<sup>\*</sup>Light Opaque Shades – Pigment content < 6wt%

\*\*Dark Opaque Shades – Pigment content > 6wt%

10

15

20

25

30

The low glass transition temperature of the opaque porcelain is paramount to assure its compatibility with alloys having CTE's in the range of about 15.5 to about 17, such as commercially available Gold Core 75<sup>TM</sup> alloy from Jeneric/Pentron Inc., Wallingford, CT. This is a gold alloy that does not contain copper and other elements that form dark oxide layers and, therefore, requires much less intricate preparation procedures compared to the Degunorm alloy. Specifically, the Gold Core 75<sup>TM</sup> alloy forms an adequate oxide layer without compromising the appearance of the coping when degassed at 870 °C-885 °C for about 5-7 min in air or vacuum.

The high potassium content in the porcelain is essential to assure high stability of leucite. High potassium oxide content combined synergetically with other alkali elements (Li and Na) assures as well relatively low glass transition temperature (GTT) and, hence, increased resistance to thermal expansion mismatch cracking and increased adaptability to alloys of slightly lower expansion. It was surprisingly found that increased potassium content increases stability in compositions with low GTT, e.g., Example 1 was found to be much more stable than Comparative Example 1 (compare K<sub>2</sub>O content). Specifically, dental porcelain of Example 1 has excellent thermal expansion stability and maintains the same thermal expansion after 5 successive bakes at its firing temperature. Dental porcelain of Comparative Example 1 was found to change thermal expansion and opacity upon multiple bakes.

Essential to this invention is that opaque compositions possess a relatively low transition temperature and contain the same elements such as ZnO and  $Ta_2O_5$  as the oxide layer forming on the alloy that assure good bonding to alloys. Specifically, the oxide layer on the Gold Core 75<sup>TM</sup> alloy was found to be enriched with ZnO and  $Ta_2O_5$  and the same components were included in the opaque formulation to improve bonding.

Body and incisal porcelain are typically applied to opaque porcelain, respectively. Preferably, the body and incisal porcelains used with the opaque herein exhibit an average coefficient of thermal expansion of about 17.2.

In a preferred embodiment of the invention, alloys having a CTE in the range of about 15.5 to about  $17 \times 10^{-6}$ °C in the temperature range of 25 ° – 500 ° C are used

to manufacture a metal core for a restoration. Opaque porcelains are applied thereto, wherein the CTE is in the range of about 16.9 to about 17.5 x  $10^{-6}$ /°C in the temperature range of 25 ° – 500 ° C and body porcelains are applied thereto having CTEs in the range of about 16.9 to about 17.7. It is preferable that the CTE of the opaque is between the CTE of the alloy and the CTE of the body porcelain.

The following examples illustrate the invention.

5

10

15

20

#### Examples

Copings and bridge frameworks cast from Bio-75G, GoldCore75, GoldCore55 and JewelCast alloys available from Jeneric/Pentron, Wallingford, CT, were prepared with a carbide tool, sand-blasted with alumina sand at pressure of 2 bar and ultrasonically cleaned in water for about 5 min. The same degassing cycle given in the firing charts below was used for Bio-75G, JewelCast, GoldCore55 and GoldCore75 castings. Following degassing, the oxide layer was removed by sand-blasting and castings were ultrasonically cleaned in water for about 5 min. The opaque of the composition of Example 3 (Table 6) was applied in two thin coats and fired according to the firing cycle given in a table below. Body/Incisal porcelain of composition of Example 1 was used to build full contour crowns and bridges and fired up to 5 times as per firing chart below.

No cracking was observed on single unit restorations made from the alloys above. However, cracks in pontic areas were found when porcelain was fired onto bridge frameworks made from Bio-75G. Both single and multiunit restorations made from GoldCore 75, GoldCore 55 and JewelCast exhibited no cracking upon multiple firings.

Firing chart in °F

	Degassing cycle	Opaque bake (2 coats)	1 <sup>st</sup> OPC Low Wear bake	2 <sup>nd</sup> OPC Low Wear bake	3 <sup>rd</sup> -5 <sup>th</sup> OPC Low Wear bake
Predry, min	0	6	6	6	6
Low T, °F	1200	600	1000	1000	1000
High T, °F	1625	1600	1575	1550	1550
Rate, °C/min	100	75	75	75	75
Vacuum	100%	100%	100%	100%	100%
VacOn, °F	1200	750	1000	1000	1000
VacOff, °F	1625	1500	1525	1500	1500
Hold, min	5 in vacuum	0	0	0	0
Cool, min	0	0	0	0	0

Firing chart in °C

	Degassing	Opaque bake	1st OPC Low	2 <sup>nd</sup> OPC Low	3 <sup>rd</sup> -5 <sup>th</sup>	OPC
	cycle	(2 coats)	Wear bake	Wear bake	Low	Wear
					bake	
Predry, min	0	6	6	6	6	
Low T, °C	650	316	538	538	538	
High T, °C	885	871	857	843	843	
Rate, °C/min	55	42	42	42	42	
Vacuum	100%	100%	100%	100%	100%	
VacOn, °C	650	399	538	538	538	
VacOff, °C	885	816	829	816	816	
Hold, min	5 in vacuum	0	0	0	0	
Cool, min	0	0	0	0	0	

In addition to dental restorations, bond flags were cast from the alloys listed above. Two thin coats of opaque (composition of Example 3) were applied and fired onto the bond flags. Bond flags were bent using pliers and metal surface exposed along the bend where opaque is fractured was inspected using optical stereomicroscope under magnification of 10x. Fracture along opaque-metal interface was found mostly adhesive, i.e. substantial fraction of the metal surface was covered by opaque indicating good bonding between alloys and opaque. The observed coverage was comparable to other metal-porcelain systems and therefore deemed sufficient. Bond strength was quantified according to ISO-9693 Metal-Ceramic Bond

Test (Schwickerath crack initiation test). The following Table indicates the bond strengths calculated from the formula  $\tau_b = k \cdot F_{fail}$ 

wherein

 $\tau_b$  is the debonding/crack initiation strength k is a coefficient which is a function of the thickness of the metal substrate, and the value of Young's modulus of the used

metallic material; and

Ffail is the fracture force

Specimen	Thick- ness (mm)	Elastic Modulus (GPa)	Load (Lbs)	F(fail) Newtons	ĸ	τb (MPa)
gold core 75	0.55	12.5	2.81	12.49950291	3.6	44.99821048
gold core 75	0.54	12.5	2.262	10.06187743	3.75	37.73204037
gold core 75	0.55	12.5	1.81	8.051281235	3.65	29.38717651
gold core 75	0.55	12.5	3.02	13.43362946	3.65	49.03274754
gold core 75	0.55	12.5	2.78	12.36605626	3.65	45.13610535
mean				0		41.25725605
Std Dev				0		7.791412161
gold core 55	0.5	15.12	1.42	6.316474781	4.1	25.8975466
gold core 55	0.5	15.12	1.07	4.759597194	4.1	19.5143485
gold core 55	0.52	15.12	2.88	12.81087843	3.8	48.68133803
gold core 55	0.5	15.12	1.68	7.473012417	4.1	30.63935091
gold core 55	0.47	15.12	1.51	6.716814732	4.7	31.56902924
gold core 55	0.47	15.12	0.98	4.359257243	4.7	20.48850904
mean		]		0		29.46502039
Std dev						9.724481694
i	1			ı		

As will be appreciated, the present invention provides porcelain compositions compatible with alloys for use in the manufacture of PFM restorations. The porcelains exhibit a coefficient of thermal expansion (CTE) substantially equal to or slightly above the CTE of the metal to which it is applied. The porcelains are fabricated from a mixture of two frit compositions. A high expansion, leucite containing frit is combined with a low melting glass frit to provide a porcelain having

10

an expansion in the range of 16.9 to about 18 x 10<sup>-6</sup>/°C in the temperature range of 25° - 500°C. By combining two frits, the expansion and firing temperature can be easily controlled.

While various descriptions of the present invention are described above, it should be understood that the various features can be used singly or in any combination thereof. Therefore, this invention is not to be limited to only the specifically preferred embodiments depicted herein.

Further, it should be understood that variations and modifications within the spirit and scope of the invention may occur to those skilled in the art to which the invention pertains. Accordingly, all expedient modifications readily attainable by one versed in the art from the disclosure set forth herein that are within the scope and spirit of the present invention are to be included as further embodiments of the present invention. The scope of the present invention is accordingly defined as set forth in the appended claims.

5

#### What Is Claimed Is:

1. A porcelain composition for use as an opaque on dental alloys in the manufacture of a dental restoration comprising by weight percent:

about 48 to about 65 % SiO<sub>2</sub>; about 10 to about 15 % Al<sub>2</sub>O<sub>3</sub>; about 0.5 to about 2 % CaO; about 1.5 to about 3 % Li<sub>2</sub>O; about 15 to about 17 % K<sub>2</sub>O; about 4 to about 6 % Na<sub>2</sub>O; and about 0.4 to about 1 F.

- 2. The porcelain of claim 1 possessing a coefficient of thermal expansion slightly higher than the dental alloys to which it is applied as the opaque.
- 3. The porcelain of claim 1 wherein the coefficient of thermal expansion is equal to or up to about  $1.5 \times 10^{-6}$ °C higher than the dental alloys to which it is applied as the opaque.
- 4. The porcelain composition of claim 1 possessing a coefficient of thermal expansion which is compatible with alloys possessing a coefficient of thermal expansion in the range of about 15.5 to about 17 x 10<sup>-6</sup>/°C in the temperature range of 25° 500°C.
- 5. The porcelain of claim 1 possessing a coefficient of thermal expansion of about 16.9 to about 18.5 in the temperature range of 25°C to 500 °C.

10

5

6. The composition of claim 1 further comprising by weight:

about 0 to about 0.7 % B<sub>2</sub>O<sub>3</sub>;
about 0 to about 5 % ZnO;
about 0 to about 2 % MgO;
about 0 to about 1 % BaO;
about 0 to about 2 % TiO<sub>2</sub>;
about 0 to about 17 % ZrO<sub>2</sub>;
about 0 to about 7 % ZrSiO<sub>4</sub>;
about 0 to about 1 % CeO<sub>2</sub>;
about 0 to about 2 % Ta<sub>2</sub>O<sub>5</sub>; and
about 0 to about 18 % SnO<sub>2</sub>.

7. A method of making a dental restoration comprising:

forming a dental porcelain powder from a dental composition comprising

about 48 to about 65 % SiO<sub>2</sub>, about 10 to about 15 % Al<sub>2</sub>O<sub>3</sub>, about 0.5 to about 2 % CaO, about 1.5 to about 3 % Li<sub>2</sub>O, about 14 to about 17 % K<sub>2</sub>O, about 4 to about 6 % Na<sub>2</sub>O, and about 0.4 to about 1 F, wherein the maturing temperature is in the range of about 750°C to about 890°C;

shaping the dental porcelain powder onto a metal core; and heating the shaped dental porcelain powder to between about 750°C to

10 about

880°C to fuse the dental porcelain powder to the metal core;

wherein the metal core exhibits a coefficient of thermal expansion in the range from about 15.5 to about  $17 \times 10^{-6}$ /°C (measured from 25°C to 500°C).

8. The method of claim 7 wherein the porcelain has a coefficient of thermal in the range of about 17 to about 17.5 in the temperature range of 25°C to 470 °C

WO 01/01925 PCT/US00/40296

15

9. The method of claim 7 wherein the metal framework comprises a gold alloy.

- 10. The method of claim 7 wherein the porcelain powder is an opaque porcelain.
- 11. The method of claim 10 further comprising applying a body porcelain over the opaque porcelain.
- 12. The method of claim 11 further comprising applying an incisal porcelain over the body porcelain.
  - 13. The method of claim 12 wherein the incisal porcelain comprises:

about 59 to about 65 % SiO<sub>2</sub>; about 10 to about 15 % Al<sub>2</sub>O<sub>3</sub>; about 0.5 to about 2 % CaO; about 1.5 to about 3 % Li<sub>2</sub>O; about 15 to about 17 % K<sub>2</sub>O; about 4 to about 6 % Na<sub>2</sub>O; and about 0.4 to about 1 F.

5

10

5

14. A porcelain composition comprising by weight percent:

about 59 to about 65 % SiO<sub>2</sub>, about 10 to about 15 % Al<sub>2</sub>O<sub>3</sub>; about 0.5 to about 2 % CaO; about 1.5 to about 3 % Li<sub>2</sub>O; about 12 to about 17 % K<sub>2</sub>O;

about 4 to about 6 % Na<sub>2</sub>O; about 0.4 to about 1 % F; and

about 10 to about 20 % of an opacifier;

wherein the porcelain composition is used as an opaque in the manufacture of dental restorations.

15. The porcelain composition of claim 16 wherein the opacifier is selected

from Al<sub>2</sub>O<sub>3</sub>, ZnO, TiO<sub>2</sub>, ZrO<sub>2</sub>, ZrSiO<sub>4</sub>, CeO<sub>2</sub>, Ta<sub>2</sub>O<sub>5</sub>, SnO<sub>2</sub> and mixtures thereof.

16. A dental restoration comprising:

a metal core, having a coefficient of thermal expansion below about 17  $\times$  10<sup>-6</sup>/°C (measured from 25°C to 500°C); and

an opaque porcelain applied on the metal core, having a coefficient of thermal expansion about equal to or up to about  $1.5 \times 10^{-6}$ /°C higher than the metal core.

5

17. The dental restoration of claim 16 wherein the opaque porcelain comprises by weight:

> about 48 to about 65 % SiO<sub>2;</sub> about 10 to about 15 % Al<sub>2</sub>O<sub>3</sub>; about 0.5 to about 2 % CaO; about 1.5 to about 3 % Li<sub>2</sub>O; about 14 to about 17 % K<sub>2</sub>O; about 4 to about 6 % Na<sub>2</sub>O; and about 0.4 to about 1 F.

18. The dental restoration of claim 17 wherein the opaque porcelain further comprises:

about 0 to about 0.7 % B<sub>2</sub>O<sub>3</sub>; about 0 to about 5 % ZnO; about 0 to about 2 % MgO; about 0 to about 1 % BaO; about 0 to about 2 % TiO2; about 0 to about 17 % ZrO<sub>2</sub>; about 0 to about 7 % ZrSiO<sub>4</sub>; 10 about 0 to about 1 % CeO<sub>2</sub>; about 0 to about 2 % Ta<sub>2</sub>O<sub>5</sub>; and about 0 to about 18 % SnO<sub>2</sub>.

10

19. A dental restoration comprising:

a metal core, having a coefficient of thermal expansion below about 17 x 10<sup>-6</sup>/°C (measured from 25°C to 500°C);

an opaque porcelain applied on the metal core, having a coefficient of thermal expansion about equal to or up to about  $1.5 \times 10^{-6}$ /°C higher than the metal core; and

a body porcelain applied to the opaque porcelain having a coefficient of thermal expansion about equal to or up to about  $1.5 \times 10^{-6}$  C higher than the metal core.

20. The dental restoration of claim 19 wherein the coefficient of thermal expansion of the opaque porcelain is between the coefficient of thermal expansion of the alloy and the coefficient of thermal expansion of the body porcelain.

# INTERNATIONAL SEARCH REPORT

Intern. Jal Application No PCT/US 00/40296

A CLASSIF	CATION OF SUBJECT MATTER A61K6/02 C03C10/00		
*	international Patent Classification (IPC) or to both national classificat	tion and IPC	
B. FIELDS S			
Minimum doc	currentation searched (classification system followed by classification A61K C03C	n aymbols)	
	on searched other than minimum documentation to the extent that su		
	ta base consulted during the international search (name of data base ternal, WPI Data, PAJ	e and, where practical, eearch terms used)	,
C. DOCUME	ENTS CONSIDERED TO BE RELEVANT		
Category *	Citation of document, with indication, where appropriate, of the rele	want peasages	Relevant to claim No.
X	US 5 653 791 A (DIMEGLIO LISA M 5 August 1997 (1997-08-05) column 2, line 9 -column 3, line column 5, line 20 - line 56 claims; tables		1,4-8, 10-13
x	EP 0 518 454 A (IVOCLAR AG) 16 December 1992 (1992-12-16) page 2, line 49 -page 3, line 2 page 3, line 17 - line 23 page 3, line 49 - line 51		1
X	EP 0 870 479 A (MEGADENTA GMBH DENTALPRODUKTE) 14 October 1998 (1998-10-14) column 4, line 10 - line 20	-/	1
X Furt	her documents are listed in the continuation of box C.	X Patent family members are listed	in annex.
"A" docum consis "E" earlier filing ( "L" docum which citatic "O" docum other "P" docum later (	ent defining the general state of the art which is not dered to be of particular relevance document but published on or after the international data ent which may throw doubts on priority claim(s) or it is cited to establish the publication date of another on or other special reason (as specified) enter referring to an oral disclosure, use, exhibition or means ent published prior to the international filing date but than the priority date claimed	To later document published after the integer or priority date and not in conflict with cited to understand the principle or the invention  "X" document of particular relevance; the common to considered novel or cannot involve an inventive step when the document of particular relevance; the common to common to common the common to the same patent.  "&" document member of the same patent.	the application but sony underlying the laimed invention be considered to purment is taken alone laimed invention rentive step when the re other auch docuse to a person sidilled family
1	lO November 2000	20/11/2000	
Name and	mailing address of the ISA  European Patent Office, P.B. 5818 Patentiasn 2  NL - 2280 HV Rijswijk  Tel. (+31-70) 340-2040, Tx. 31 651 epo nl, Fax: (+31-70) 340-3016	Authorized officer  Cousins-Van Steen	, G

## INTERNATIONAL SEARCH REPORT

Intern. Jail Application No PCT/US 00/40296

		PC1/05 00	7 40230
C.(Continue	INTERPOLATION DOCUMENTS CONSIDERED TO BE RELEVANT		
Category *	Citation of document, with indication, where appropriate, of the relevant passages		Relevant to claim No.
A	WO 99 18910 A (JENERIC PENTRON INC; PANZERA CARLINO (US)) 22 April 1999 (1999-04-22) page 5, line 11 -page 7, line 3 page 11, line 1 - line 10		1-13
A	WO 95 11866 A (MINNESOTA MINING & MFG) 4 May 1995 (1995-05-04)		
A	US 5 453 290 A (VAN DER ZEL JOSEPH M) 26 September 1995 (1995-09-26) cited in the application		
A	EP 0 272 745 A (AMERICAN THERMOCRAFT CORP) 29 June 1988 (1988-06-29)		

# INTERNATIONAL SEARCH REPORT

information on patent family members

Intern. hel Application No PCT/US 00/40296

Patent document cited in search report		Publication date	Patent family member(s)		Publication date
US 5653791	Α	05-08-1997	AU	1621397 A	22-01-1998
00 0000/51	••		CA	2199569 A	12-09-1997
			EP	0795311 A	17-09-1997
			JP	10036137 A	10-02-1998
			US	5944884 A	31-08-1999
EP 0518454	A	16-12-1992	DE	4119483 A	24-12-1992
	•	<b>30 12 3</b>	AT	126049 T	15-08-1999
			AU	650338 B	16-06-1994
			AU	1739592 A	17-12-199
			CA	2070691 A	14-12-1992
			DE	59203186 D	14-09-199
EP 0870479	A	14-10-1998	DE	19714178 A	08-10-199
W0 9918910	A	22-04-1999	AU	9802698 A	03-05-199
WO 9511866	A	04-05-1995	US	5622551 A	22-04-199
US 5453290	A	26-09-1995	NL	9001986 A	01-04-199
	••	20 00 2000	AT	110254 T	15-09-199
			DE	475528 T	25-02-199
			DK	475528 T	19-09-199
			EP	0475528 A	18-03-199
			ES	2059042 T	01-11-199
			JP	4230611 A	19-08-199
EP 0272745	Α	29-06-1988	US	4798536 A	17-01-198
	••		AT	75137 T	15-05-199
			CA	1272222 A	31-07-199
			DE	3778513 A	27-05-199
			JP	63287709 A	24-11-198
			JP	1882859 C	10-11-199
			JP	6000662 B	05-01-199
			JP	63265856 A	02-11-198